

**Copper base, Type: Bronze - Phosphorus, RM, Mushroom 60 mm x 5 mm**

Art. ID                      CTIF-B12  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,12	%		
	Arsenic (As)	[7440-38-2]	0,111	%		
	Copper (Cu)	[7440-50-8]	85,65	%		
	Iron (Fe)	[7439-89-6]	0,162	%		
	Manganese (Mn)	[7439-96-5]	0,235	%		
	Nickel (Ni)	[7440-02-0]	2,63	%		
	Phosphorus (P)	[7723-14-0]	0,525	%		
	Lead (Pb)	[7439-92-1]	0,201	%		
	Sulfur (S)	[7704-34-9]	0,013	%		
	Antimony (Sb)	[7440-36-0]	0,117	%		
	Silicon (Si)	[7440-21-3]	0,05	%		
	Tin (Sn)	[7440-31-5]	9,57	%		
	Zinc (Zn)	[7440-66-6]	0,61	%		