

Copper base, Type: Bronze - Tin, RM, 936, Mushroom 60 mm x 5 mm

Art. ID CTIF-B31
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	78,65	%		
	Nickel (Ni)	[7440-02-0]	0,489	%		
	Lead (Pb)	[7439-92-1]	11,79	%		
	Sulfur (S)	[7704-34-9]	0,028	%		
	Antimony (Sb)	[7440-36-0]	0,475	%		
	Tin (Sn)	[7440-31-5]	7,65	%		
	Zinc (Zn)	[7440-66-6]	0,79	%		