

Copper base, Type: Bronze - Tin, RM, 945 MOD, Mushroom 60 mm x 5 mm

Art. ID CTIF-B32
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,075	%		
	Arsenic (As)	[7440-38-2]	0,0056	%		
	Copper (Cu)	[7440-50-8]	74,8	%		
	Iron (Fe)	[7439-89-6]	0,11	%		
	Nickel (Ni)	[7440-02-0]	1,49	%		
	Phosphorus (P)	[7723-14-0]	0,039	%		
	Lead (Pb)	[7439-92-1]	16,1	%		
	Sulfur (S)	[7704-34-9]	0,027	%		
	Antimony (Sb)	[7440-36-0]	0,13	%		
	Silicon (Si)	[7440-21-3]	0,07	%		
	Tin (Sn)	[7440-31-5]	5,92	%		
	Zinc (Zn)	[7440-66-6]	1,17	%		