

Copper base, Type: Bronze - Leaded, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-B4
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	83,75	%		
	Iron (Fe)	[7439-89-6]	0,021	%		
	Nickel (Ni)	[7440-02-0]	0,57	%		
	Phosphorus (P)	[7723-14-0]	0,52	%		
	Lead (Pb)	[7439-92-1]	2,53	%		
	Sulfur (S)	[7704-34-9]	0,019	%		
	Antimony (Sb)	[7440-36-0]	0,1	%		
	Silicon (Si)	[7440-21-3]	0,015	%		
	Tin (Sn)	[7440-31-5]	11,1	%		
	Zinc (Zn)	[7440-66-6]	1,34	%		