

**Anti-ORAI2 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA065937-25UL  
Unit                         1 x 25 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-C7orf19|Anti-CBCIP2|Anti-FLJ12474|Anti-FLJ14733|Anti-H\_NH0514P08.8|Anti-TMEM142B|