

Copper base, Type: Bronze - Aluminum, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-CA37
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	9,84	%		
	Bismuth (Bi)	[7440-69-9]	0,0118	%		
	Chromium (Cr)	[7440-47-3]	0,085	%		
	Copper (Cu)	[7440-50-8]	76,79	%		
	Iron (Fe)	[7439-89-6]	6,85	%		
	Magnesium (Mg)	[7439-95-4]	0,077	%		
	Manganese (Mn)	[7439-96-5]	0,752	%		
	Nickel (Ni)	[7440-02-0]	4,98	%		
	Lead (Pb)	[7439-92-1]	0,0503	%		
	Silicon (Si)	[7440-21-3]	0,04	%		
	Tin (Sn)	[7440-31-5]	0,147	%		
	Zinc (Zn)	[7440-66-6]	0,364	%		