

Copper base, Type: Nickel Alloy, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-CN33
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	0,0212	%		
	Carbon (C)	[7440-44-0]	0,02	%		
	Cadmium (Cd)	[7440-43-9]	0,006	%		
	Copper (Cu)	[7440-50-8]	66,9	%		
	Iron (Fe)	[7439-89-6]	1,6	%		
	Magnesium (Mg)	[7439-95-4]	<0,06	%		
	Manganese (Mn)	[7439-96-5]	0,45	%		
	Niobium (Nb)	[7440-03-1]	0,06	%		
	Nickel (Ni)	[7440-02-0]	29,75	%		
	Phosphorus (P)	[7723-14-0]	0,02	%		
	Lead (Pb)	[7439-92-1]	0,053	%		
	Sulfur (S)	[7704-34-9]	0,013	%		
	Silicon (Si)	[7440-21-3]	0,47	%		
	Tellurium (Te)	[13494-80-9]	0,0224	%		
	Zinc (Zn)	[7440-66-6]	0,37	%		