

Copper base, Type: Nickel Alloy, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-CN4
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	84	%		
	Iron (Fe)	[7439-89-6]	1,8	%		
	Manganese (Mn)	[7439-96-5]	1,5	%		
	Niobium (Nb)	[7440-03-1]	0,7	%		
	Nickel (Ni)	[7440-02-0]	11,2	%		
	Lead (Pb)	[7439-92-1]	0,006	%		
	Sulfur (S)	[7704-34-9]	~0,001	%		
	Tin (Sn)	[7440-31-5]	0,058	%		
	Zinc (Zn)	[7440-66-6]	0,07	%		