

Copper base, Type: Beryllium-Cobalt, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-CuBeCo6
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1,37	%		
	Aluminium (Al)	[7429-90-5]	0,135	%		
	Beryllium (Be)	[7440-41-7]	1,54	%		
	Cobalt (Co)	[7440-48-4]	<2	%		
	Chromium (Cr)	[7440-47-3]	0,0576	%		
	Copper (Cu)	[7440-50-8]	93,09	%		
	Iron (Fe)	[7439-89-6]	0,109	%		
	Manganese (Mn)	[7439-96-5]	0,0173	%		
	Nickel (Ni)	[7440-02-0]	<1,5	%		
	Lead (Pb)	[7439-92-1]	0,0397	%		
	Silicon (Si)	[7440-21-3]	0,26	%		
	Tin (Sn)	[7440-31-5]	0,0135	%		
	Zinc (Zn)	[7440-66-6]	0,033	%		