

Anti-ATP6V0A4 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA064555-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms:

|Anti-ATP6N1B|Anti-ATP6N2|Anti-RDRTA2|Anti-RTA1C|Anti-RTADR|Anti-Stv1|Anti-VPP2|Anti-Vph1|Anti-a4|