

**Copper base, Type: Brass, RM, Mushroom 60 mm x 5 mm**

Art. ID CTIF-L21  
Unit each  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,103	%		
	Copper (Cu)	[7440-50-8]	82,5	%		
	Iron (Fe)	[7439-89-6]	0,086	%		
	Manganese (Mn)	[7439-96-5]	0,004	%		
	Nickel (Ni)	[7440-02-0]	0,156	%		
	Phosphorus (P)	[7723-14-0]	0,05	%		
	Lead (Pb)	[7439-92-1]	0,209	%		
	Antimony (Sb)	[7440-36-0]	0,1	%		
	Silicon (Si)	[7440-21-3]	0,036	%		
	Tin (Sn)	[7440-31-5]	1,5	%		
	Zinc (Zn)	[7440-66-6]	15,4	%		