

Copper base, Type: Brass, RM, Mushroom 60 mm x 5 mm

Art. ID CTIF-L23
Unit each
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,051	%		
	Copper (Cu)	[7440-50-8]	81,2	%		
	Iron (Fe)	[7439-89-6]	0,246	%		
	Nickel (Ni)	[7440-02-0]	0,033	%		
	Lead (Pb)	[7439-92-1]	0,058	%		
	Silicon (Si)	[7440-21-3]	0,28	%		
	Tin (Sn)	[7440-31-5]	0,2	%		
	Zinc (Zn)	[7440-66-6]	17,9	%		