

Anti-IL37 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID SAF-HPA057950-25UL
Unit 1 x 25 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-FIL1|Anti-FIL1(ZETA)|Anti-FIL1Z|Anti-IL-1F7|Anti-IL-1H4|Anti-IL1F7|Anti-IL1F7|