

**Copper base, Type: Brass - Silicon, RM, Mushroom 60 mm x 5 mm**

Art. ID                      CTIF-LS1  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,02	%		
	Copper (Cu)	[7440-50-8]	77,7	%		
	Iron (Fe)	[7439-89-6]	0,448	%		
	Manganese (Mn)	[7439-96-5]	0,039	%		
	Nickel (Ni)	[7440-02-0]	0,55	%		
	Phosphorus (P)	[7723-14-0]	0,128	%		
	Lead (Pb)	[7439-92-1]	0,213	%		
	Silicon (Si)	[7440-21-3]	4,35	%		
	Tin (Sn)	[7440-31-5]	0,243	%		
	Zinc (Zn)	[7440-66-6]	16,3	%		