

Anti-PTP4A1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID	SAF-HPA003281-100UL
Unit	1 x 100 µL
Deliverydetails	/not restricted

Description

Synonyms: |Anti-PRL-1|Anti-PTP(CAAXI)|Anti-Protein tyrosine phosphatase type IVA 1|Anti-Protein- tyrosine phosphatase 4a1|Anti-Protein-tyrosine phosphatase of regenerating liver 1|