

Anti-SAFB antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution

Art. ID	SAF-HPA020076-100UL
Unit	1 x 100 µL
Deliverydetails	/not restricted

Description

Synonyms: |Anti-HSP27 ERE-TATA-binding protein|Anti-HSP27 estrogen response element-TATA box-binding protein|Anti-SAB-B1|Anti-SAF-B|Anti-Scaffold attachment factor B1|