

**Anti-POGLUT1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution**

Art. ID                      SAF-HPA037855-100UL  
Unit                         1 x 100 µL  
Deliverydetails            /not restricted

Description

Synonyms: |Anti-9630046K23Ri|Anti-C3orf9|Anti-HCLP46|Anti-KDELCL1|Anti-KTEL (Lys-Tyr-Glu-Leu) containing 1|Anti-KTELC1|Anti-MDS010|Anti-MDSRP|Anti-MGC32995|