

Anti-KYAT1 antibody produced in rabbit, Prestige Antibodies(R) Powered by Atlas Antibodies, affinity isolated antibody, buffered aqueous glycerol solution, Ab1

Art. ID SAF-HPA021176-100UL
Unit 1 x 100 µL
Deliverydetails /not restricted

Description

Synonyms: |Anti-CCBL1|Anti-GTK|Anti-KATI|Anti-Cysteine-S-conjugate
beta-lyase|Anti-GTK|Anti-Glutamine--phenylpyruvate transaminase|Anti-KATI|Anti-Kynurenine aminotransferase
I|Anti-Kynurenine--oxoglutarate transaminase 1|