

**Digitoxigenin-bisdigitoxoside (Chemical purity: 98.5 % (NMR), 98.1 % (HPLC))**

Art. ID LM24-AC-10014  
Unit on request  
Deliverydetails De Minimis

| Text/Information | Analyte/Parameter                 | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------------------|--------------|---------------------|------|--------|--------|
|                  | Digitoxigenin-bisdigito<br>xoside | [16479-50-8] |                     |      |        |        |