

Custom Solution - 8 components in HNO3 7%

Art. ID LM24-CUS-71179
Unit 2 x 500 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Description

Conservation 12 months. (6 months after opening)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 10 | mg/L | | |
| | Sulfur (S) | [7704-34-9] | 2000 | mg/L | | |
| | Potassium (K) | [7440-09-7] | 500 | mg/L | | |
| | Silicon (Si) | [7440-21-3] | 250 | mg/L | | |
| | Magnesium (Mg) | [7439-95-4] | 500 | mg/L | | |
| | Calcium (Ca) | [7440-70-2] | 500 | mg/L | | |
| | Sodium (Na) | [7440-23-5] | 1000 | mg/L | | |
| | Phosphorus (P) | [7723-14-0] | 100 | mg/L | | |