

Naval Brass, CDA 486 / C48600, Ø 38x19mm

Art. ID IARM-CU486-18-G
Unit ea

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Certified Value	Silver (Ag)	[7440-22-4]	0,004 ± 0,002	%		
Certified Value	Arsenic (As)	[7440-38-2]	0,025 ± 0,007	%		
Certified Value	Cadmium (Cd)	[7440-43-9]	0,0009 ± 0,0002	%		
Certified Value	Copper (Cu)	[7440-50-8]	61,2 ± 0,3	%		
Certified Value	Iron (Fe)	[7439-89-6]	0,036 ± 0,003	%		
Certified Value	Nickel (Ni)	[7440-02-0]	0,032 ± 0,003	%		
Certified Value	Phosphorus (P)	[7723-14-0]	0,004 ± 0,002	%		
Certified Value	Lead (Pb)	[7439-92-1]	1,31 ± 0,05	%		
Certified Value	Tin (Sn)	[7440-31-5]	0,692 ± 0,008	%		
Certified Value	Zinc (Zn)	[7440-66-6]	36,5 ± 0,3	%		
Indicative Value	Aluminium (Al)	[7429-90-5]	~30	ppm		
Indicative Value	Boron (B)	[7440-42-8]	<50	ppm		
Indicative Value	Bismuth (Bi)	[7440-69-9]	~4	ppm		
Indicative Value	Carbon (C)	[7440-44-0]	~16	ppm		
Indicative Value	Cobalt (Co)	[7440-48-4]	~6	ppm		
Indicative Value	Chromium (Cr)	[7440-47-3]	<50	ppm		
Indicative Value	Hydrogen (H)	[1333-74-0]	<10	ppm		
Indicative Value	Hafnium (Hf)	[7440-58-6]	<10	ppm		
Indicative Value	Magnesium (Mg)	[7439-95-4]	<50	ppm		
Indicative Value	Manganese (Mn)	[7439-96-5]	~3	ppm		
Indicative Value	Molybdenum (Mo)	[7439-98-7]	<50	ppm		
Indicative Value	Nitrogen (N)	[7727-37-9]	<10	ppm		
Indicative Value	Niobium (Nb)	[7440-03-1]	<50	ppm		
Indicative Value	Oxygen (O)	[7782-44-7]	~7	ppm		
Indicative Value	Sulfur (S)	[7704-34-9]	~30	ppm		
Indicative Value	Antimony (Sb)	[7440-36-0]	~50	ppm		
Indicative Value	Selenium (Se)	[7782-49-2]	~40	ppm		
Indicative Value	Silicon (Si)	[7440-21-3]	~20	ppm		
Indicative Value	Titanium (Ti)	[7440-32-6]	~10	ppm		
Indicative Value	Vanadium (V)	[7440-62-2]	<50	ppm		
Indicative Value	Tungsten (W)	[7440-33-7]	<10	ppm		
Indicative Value	Zirconium (Zr)	[7440-67-7]	<50	ppm		

