

### **Custom Solution - 9 components 100 mg/L each of in HNO3 5%**

Art. ID LM24-CUS-71194  
Unit 2 x 100 mL bottle  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 II E2

#### **Description**

Conservation 24 months (originally closed bottle), (12 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Bismuth (Bi)	[7440-69-9]	100	mg/L		
	Tellurium (Te)	[13494-80-9]	100	mg/L		
	Niobium (Nb)	[7440-03-1]	100	mg/L		
	Gallium (Ga)	[7440-55-3]	100	mg/L		
	Tantalum (Ta)	[7440-25-7]	100	mg/L		
	Germanium (Ge)	[7440-56-4]	100	mg/L		
	Tin (Sn)	[7440-31-5]	100	mg/L		
	Antimony (Sb)	[7440-36-0]	100	mg/L		
	Tungsten (W)	[7440-33-7]	100	mg/L		