

## **Aroclor 1254 - 50 mg/kg in Transformatoröl**

Art. ID LM24-CUS-10207  
Unit 1 x 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

### Description

Haltbarkeit 18 Monate (ungeöffnet)

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|--------------|---------------------|-------|--------|--------|
|                  | Aroclor 1254      | [11097-69-1] | 50                  | mg/kg |        |        |