

## **Aroclor 1260 - 50 µg/mL in Transformatoröl**

Art. ID LM24-CUS-10208  
Unit 1 x 5 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

### Description

Haltbarkeit 18 Monate (ungeöffnet)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1260	[11096-82-5]	50	µg/mL		