

**L-PHENYLALANINE (4'-13C, 99%, 2,3,3,2',3',5',6'-D7, 98%, 15N, 98%)**

Art. ID CIL-CDNLM-11149-PK  
Unit on request

| Text/Information | Analyte/Parameter                                                  | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------------------------------------------|------------|---------------------|------|--------|--------|
|                  | L-PHENYLALANINE (4'-13C, 99%, 2,3,3,2',3',5',6'-D7, 98%, 15N, 98%) |            |                     |      |        |        |