

**REVERSE 3,3',5-TRIIODO-L-THYRONINE:HCL (REV T3) (DIIODOPHENYL-RING-13C6,99%)
(95% CHEMICAL PURITY)**

Art. ID CIL-CLM-10601-0.005

Unit 5 mg

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	REVERSE 3,3',5-TRIIODO- L-THYRONINE:HCL (REV) (DIIODOPHENYL-RING-1 C6,99%)	[1217676-14-6]				