

**REVERSE 3,3',5-TRIIODO-L-THYRONINE:HCL (REV T3) (DIIODOPHENYL-RING-13C6,99%),  
100 µg/mL 0.1N NH3/MEOH**

Art. ID CIL-CLM-10601-C

Unit 1 mL

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	REVERSE 3,3',5-TRIIODO- L-THYRONINE:HCL (REV ) (DIIODOPHENYL-RING-1 C6,99%)	[1217676-14-6]	100	µg/mL		