

**L-PHENYLALANINE-N-T-BOC (1-13C, 99%)**

Art. ID CIL-CLM-2170-0.5  
Unit 0,5 g

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                     | CAS number   | Concentration/Value | Unit | Method | Source |
|------------------|---------------------------------------|--------------|---------------------|------|--------|--------|
|                  | L-PHENYLALANINE-N-T-B<br>(1-13C, 99%) | [84771-22-2] |                     |      |        |        |