

**POTASSIUM PHOSPHOENOL PYRUVATE (2-13C, 99%)**

Art. ID CIL-CLM-3551-PK  
Unit on request

Description

Chemical purity: 97%

| Text/Information | Analyte/Parameter                           | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|---|---------------|---------------------|------|--------|--------|
|                  | POTASSIUM PHOSPHOEI<br>YRUVATE (2-13C, 99%) | [201996-53-4] |                     |      |        |        |