

**L-3,3'-DIODOTHYRONINE (T2) (PHENOXY-13C6, 99%) 100 µg/mL IN 0.1N AMMONIA IN METHANOL**

Art. ID CIL-CLM-9222-C  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Chemical purity: 97%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-3,3'-DIODOTHYRONINE [1217459-13-6] (T2) (PHENOXY-13C6, 99%) )		100	µg/mL		