

**DL-3-HYDROXYKYNURENINE (CHEMICAL PURITY 95%) (1,2,3-13C3, 98%,  
ALPHA-AMINO-15N, 98%)**

Art. ID CIL-CNLM-10399-0.001

Unit 1 mg

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-3-HYDROXYKYNUREN [2267306-16-9] (1,2,3-13C3, 98%, ALPHA -AMINO-15N, 98%)					