

**DL-3-HYDROXYKYNURENINE (CHEMICAL PURITY 95%) (1,2,3-13C3, 98%,
ALPHA-AMINO-15N, 98%)**

Art. ID CIL-CNLM-10399-PK

Unit on request

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	DL-3-HYDROXYKYNUREN [2267306-16-9]					
	(1,2,3-13C3, 98%, ALPHA					
	-AMINO-15N, 98%)					