

**BETA-ALANINE (3-13C, 99%, 15N, 98%)**

Art. ID CIL-CNLM-3440-PK  
Unit on request

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter        | CAS number    | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------|---------------|---------------------|------|--------|--------|
|                  | BETA-ALANINE (3-13C, 99% | [404839-27-6] |                     |      |        |        |
|                  | %, 15N, 98%)             |               |                     |      |        |        |