

L-HOMOARGININE:HCL (13C7, 98%+, 15N4, 98%+)

Art. ID CIL-CNLM-4645-0.01
Unit 10 mg

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------------|------------|---------------------|------|--------|--------|
| | L-HOMOARGININE:HCL (1 [2483830-23-3] | | | | | |
| | 7, 98%+, 15N4, 98%+) | | | | | |