

1.5 mM HIS6-GB1+10 mM PBS+0.02% NaN3+0.1 mM DSS 90%(V/V) H2O/D2O (U-13C, 99%, U-15N, 99%)

Art. ID CIL-CNLM-8663-1

Unit 1,5 MM TUBE

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1.5 mM HIS6-GB1+10 mM BS+0.02% NaN3+0.1 mM I S 90%(V/V) H2O/D2O (U-1 3C, 99%, U-15N, 99%)					