

**Copper base, Type: Copper, CRM, Pins 50g of 0.31 g pins**

Art. ID IMN-001  
Unit 50 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,00271	%		
	Arsenic (As)	[7440-38-2]	0,000313	%		
	Bismuth (Bi)	[7440-69-9]	0,000139	%		
	Cadmium (Cd)	[7440-43-9]	~0,000011	%		
	Cobalt (Co)	[7440-48-4]	0,000209	%		
	Chromium (Cr)	[7440-47-3]	0,000353	%		
	Iron (Fe)	[7439-89-6]	0,00139	%		
	Manganese (Mn)	[7439-96-5]	0,000516	%		
	Nickel (Ni)	[7440-02-0]	0,000364	%		
	Lead (Pb)	[7439-92-1]	0,000402	%		
	Sulfur (S)	[7704-34-9]	0,000599	%		
	Antimony (Sb)	[7440-36-0]	0,00024	%		
	Selenium (Se)	[7782-49-2]	0,000134	%		
	Tin (Sn)	[7440-31-5]	0,000339	%		
	Tellurium (Te)	[13494-80-9]	0,000123	%		
	Zinc (Zn)	[7440-66-6]	0,000457	%		