

**Copper base, Type: Copper, CRM, Refined Copper, Chips 200 g**

Art. ID IMN-4  
Unit 200 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,21	%		
	Arsenic (As)	[7440-38-2]	0,054	%		
	Bismuth (Bi)	[7440-69-9]	~0,0011	%		
	Cobalt (Co)	[7440-48-4]	0,001	%		
	Chromium (Cr)	[7440-47-3]	<0,00005	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,001	%		
	Manganese (Mn)	[7439-96-5]	<0,00003	%		
	Nickel (Ni)	[7440-02-0]	0,078	%		
	Lead (Pb)	[7439-92-1]	0,2	%		
	Antimony (Sb)	[7440-36-0]	~0,0041	%		
	Silicon (Si)	[7440-21-3]	~0,00029	%		
	Tin (Sn)	[7440-31-5]	0,00044	%		
	Zinc (Zn)	[7440-66-6]	~0,0017	%		