

**PHOSPHORYLCHOLINE CHLORIDE:4H2O CA SALT (TRIMETHYL-D9, 98%) (90%  
CHEMICAL PURITY)**

Art. ID CIL-DLM-298-PK

Unit on request

| Text/Information | Analyte/Parameter      | CAS number     | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|----------------|---------------------|------|--------|--------|
|                  | PHOSPHORYLCHOLINE C    | [2483831-33-8] |                     |      |        |        |
|                  | IDE:4H2O CA SALT (TRIM |                |                     |      |        |        |
|                  | THYL-D9, 98%)          |                |                     |      |        |        |