

**L-ARGININE:HCL (<5% D) (4,4,5,5-D4, 94%)**

Art. ID CIL-DLM-6038-PK  
Unit on request

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter      | CAS number     | Concentration/Value | Unit | Method | Source |
|------------------|------------------------|----------------|---------------------|------|--------|--------|
|                  | L-ARGININE:HCL (<5% D) | [2483831-55-4] |                     |      |        |        |
|                  | (4,4,5,5-D4, 94%)      |                |                     |      |        |        |