

SODIUM ESTRONE 3-SULFATE (2,4,16,16-D4, 98%) STABILIZED WITH 50% W/W TRIS

Art. ID CIL-DLM-7456-PK
Unit on request

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-----------------------|---------------|---------------------|------|--------|--------|
| | SODIUM ESTRONE 3-SUL | [285979-80-8] | | | | |
| | E (2,4,16,16-D4, 98%) | | | | | |