

**N-ISOPROPYLACRYLAMIDE (ISOPROPYL-D7, 98%)**

Art. ID CIL-DLM-7630-0.5  
Unit 0,5 g

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                    | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--------------------------------------|------------|---------------------|------|--------|--------|
|                  | N-ISOPROPYLACRYLAMIDE [1219803-32-3] |            |                     |      |        |        |
|                  | ISOPROPYL-D7, 98%)                   |            |                     |      |        |        |