

**(RS)-2-HYDROXYGLUTARIC ACID, DISODIUM SALT (2,3,3-D3, OD, 98%) (CHEMICAL PURITY 95%)**

Art. ID CIL-DLM-9104-PK  
Unit on request

| Text/Information | Analyte/Parameter                                              | CAS number     | Concentration/Value | Unit | Method | Source |
|------------------|----------------------------------------------------------------|----------------|---------------------|------|--------|--------|
|                  | (RS)-2-HYDROXYGLUTARIC ACID, DISODIUM SALT (2,3,3-D3, OD, 98%) | [2483831-91-8] |                     |      |        |        |