

BILE ACID STANDARD 2 MIX (CONJUGATED) (UNLABELLED) STANDARD

Art. ID	CIL-MSK-BA2-US-1
Unit	1 VIAL
Deliverydetails	No Dangerous Good /not restricted

Description

Reconstituting the dried down mix in 1 mL of solvent will yield an equimolar concentration of approximately 100 µM

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Glycochenodeoxycholic Acid Sodium Salt	[16564-43-5]	47	µg/vial		
	Glycocholic acid hydrate	[1192657-83-2]	48	µg/vial		
	Sodium glycodeoxycholate	[16409-34-0]	47	µg/vial		
	Lithocholylglycine	[474-74-8]	43	µg/vial		
	Glycoursodeoxycholic Acid	[64480-66-6]	45	µg/vial		
	Taurochenodeoxycholic Acid Sodium Salt	[6009-98-9]	52	µg/vial		
	Taurocholic acid sodium salt	[145-42-6]	54	µg/vial		
	Sodium Taurodeoxycholate Hydrate	[207737-97-1]	54	µg/vial		
	Taurolithocholic Acid Sodium Salt	[6042-32-6]	51	µg/vial		
	Tauroursodeoxycholic acid dihydrate	[117609-50-4]	54	µg/vial		