

**L-ARGININE:HCL (GUANIDO-15N2, 98%+)**

Art. ID CIL-NLM-395-1  
Unit 1 g

Description

Chemical purity: 98%

| Text/Information | Analyte/Parameter                                    | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|------------------------------------------------------|------------|---------------------|------|--------|--------|
|                  | L-ARGININE:HCL (GUANII [204633-92-1]<br>-15N2, 98%+) |            |                     |      |        |        |