

**3,3',5-TRIIODO-L-THYRONINE:HCL (T3) (UNLABELLED) 100 µg/mL IN 0.1N NH3 IN METHANOL (CHEMICAL PURITY 95%)**

Art. ID CIL-ULM-10573-C  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3,3',5-TRIIODO-L-THYRONINE:HCL (T3)		100	µg/mL		