

REVERSE 3,3',5-TRIIODO-L-THYRONINE:HCL (REV T3) (UNLABELLED) 100 µg/mL IN 0.1N NH3 IN METHANOL (95% CHEMICAL PURITY)

Art. ID CIL-ULM-10602-C
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Chemical purity: 95%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	REVERSE 3,3',5-TRIIODO-L-THYRONINE:HCL (REV T3) (UNLABELLED)		100	µg/mL		