

**Copper base, Type: Bronze - Leaded, CRM, Set IMN-BI1 - IMN-BI4 only, see IMN-BI1-BI4, Disc 40 mm x 25 mm**

Art. ID                      IMN-BI2  
Unit                         disc  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,077	%		
	Arsenic (As)	[7440-38-2]	0,11	%		
	Bismuth (Bi)	[7440-69-9]	0,07	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,31	%		
	Manganese (Mn)	[7439-96-5]	0,15	%		
	Nickel (Ni)	[7440-02-0]	1,46	%		
	Phosphorus (P)	[7723-14-0]	0,59	%		
	Lead (Pb)	[7439-92-1]	5,39	%		
	Sulfur (S)	[7704-34-9]	~0,0055	%		
	Antimony (Sb)	[7440-36-0]	0,43	%		
	Silicon (Si)	[7440-21-3]	0,13	%		
	Tin (Sn)	[7440-31-5]	4,18	%		
	Zinc (Zn)	[7440-66-6]	5,73	%		