

**L-3,3'-DIODOTHYRONINE (T2) 100 µg/mL IN 0.1N AMMONIA IN METHANOL
(UNLABELLED)**

Art. ID CIL-ULM-9223-C
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	L-3,3'-DIODOTHYRONINE (T2)		100	µg/mL		