

High Nitrogen Stainless Steel, 50 g Chips

Art. ID MBH-13X42027A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	0,294 ± 0,006	%		
	Silicon (Si)	[7440-21-3]	0,544 ± 0,006	%		
	Sulfur (S)	[7704-34-9]	0,0005 ± 0,0001	%		
	Phosphorus (P)	[7723-14-0]	0,0139 ± 0,0005	%		
	Manganese (Mn)	[7439-96-5]	0,356 ± 0,005	%		
	Nickel (Ni)	[7440-02-0]	0,163 ± 0,004	%		
	Chromium (Cr)	[7440-47-3]	15,25 ± 0,05	%		
	Molybdenum (Mo)	[7439-98-7]	0,99 ± 0,008	%		
	Copper (Cu)	[7440-50-8]	0,0357 ± 0,0013	%		
	Cobalt (Co)	[7440-48-4]	0,0191 ± 0,0009	%		
	Vanadium (V)	[7440-62-2]	0,048 ± 0,002	%		
	Niobium (Nb)	[7440-03-1]	0,004 ± 0,001	%		
	Tungsten (W)	[7440-33-7]	0,0192 ± 0,0012	%		
	Aluminium (Al)	[7429-90-5]	0,004 ± 0,001	%		
	Titanium (Ti)	[7440-32-6]	~0,002	%		
	Tin (Sn)	[7440-31-5]	0,0026 ± 0,0003	%		
	Nitrogen (N)	[7727-37-9]	0,402 ± 0,007	%		