

High-Ni Calibration series, 50 g Chips

Art. ID MBH-14XFeNi35D-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Provisional value	Carbon (C)	[7440-44-0]	0,035	%		
Provisional value	Silicon (Si)	[7440-21-3]	0,25	%		
Provisional value	Sulfur (S)	[7704-34-9]	0,145	%		
Provisional value	Phosphorus (P)	[7723-14-0]	0,04	%		
Provisional value	Manganese (Mn)	[7439-96-5]	0,3	%		
Provisional value	Nickel (Ni)	[7440-02-0]	34,2	%		
Provisional value	Chromium (Cr)	[7440-47-3]	0,4	%		
Provisional value	Cobalt (Co)	[7440-48-4]	0,39	%		
Provisional value	Copper (Cu)	[7440-50-8]	0,037	%		
Provisional value	Aluminium (Al)	[7429-90-5]	~0,003	%		